

# STEREO OPTICAL PRESCREENING RECORD FORM FOR THE OPTEC<sup>®</sup> VISION TESTER

NAME: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_ AGE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_ WEARING GLASSES: Y N  
 ADMINISTRATOR: \_\_\_\_\_ CONTACT LENSES: Y N

**PERIPHERAL TEST** L = 85° 70° 55° 45° NASAL 45° R = 85° 70° 55°

**DIAL AT 01: ACUITY - FAR** Score: 20/\_\_\_\_\_

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	ZN	RO	HK
2 20/100	RKS	HNC	ZOD
3 20/70	HCDV	SKZO	RNDS
4 20/50	ZROD	NSCH	VZKN
5 20/40	KHSC	OZNR	DNVC
6 20/30	ONRZV	DKHCS	KDSON
7 20/20	SDCHN	VRZKO	HSNRD

**DIAL AT 03: ACUITY - NEAR** Score: 20/\_\_\_\_\_

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	SVC	NRK	HZO
2 20/70	RNZH	DOKV	CSZN
3 20/50	CKVD	SNZR	DOHC
4 20/40	VHRN	ODSK	NZCS
5 20/30	HSKRC	NZDOV	ZSHNK
6 20/20	ZONVR	HCSKD	VKCD S

**DIAL AT 02: COLOR PERCEPTION - FAR**    A-12    B-5    C-26    D-6    E-16    F-BLANK    PASS \_\_\_\_ FAIL \_\_\_\_

**DIAL AT 04: LATERAL PHORIA - FAR**    1 2 3    4 5 6    7 8 9    10 11 12    13 14 15

1 PRISM DIOPTERS    ESOPHORIA    ORTHO    EXOPHORIA

**DIAL AT 05: VERTICAL PHORIA - FAR**    1 2    3 4 5    6 7

1/2 PRISM DIOPTERS    LEFT HYPERPHORIA    ORTHO    RIGHT HYPERPHORIA

**DIAL AT 06: LATERAL PHORIA - NEAR**    1 2 3    4 5 6    7 8 9    10 11 12    13 14 15

1 PRISM DIOPTERS    ESOPHORIA    ORTHO    EXOPHORIA

**DIAL AT 07: DEPTH PERCEPTION - FAR**    1-B    2-L    3-B    4-T    5-T    6-L    7-R    8-L    9-R

SECONDS OF ARC    400    200    100    70    50    40    30    25    20

**DIAL AT 08: FUSION - FAR**    3 cubes BOTH eyes = **PASS**    Other than 3 cubes = **FAIL**

<b>DIAL AT 09: LEFT EYE</b>		Score: 20/_____ Right Occluder OFF								<b>RIGHT EYE</b>		Score: 20/_____ Left Occluder OFF							
<b>TUMBLING</b>		<u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u>										<u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u> <u>20</u>							
<b>"E" FAR</b>		PILOT    200   100   70   50   40   30   20										PILOT    200   100   70   50   40   30   20							
<b>LINE</b>		1   2   3   4   5   6										1   2   3   4   5   6							
LINE/ACUITY	COLOR	L	T	B	L	R	L	T			R	B	R	L	L	T	B		
A 20/50	BLUE				R	T	R	L	T	B	A 20/50			L	R	T	L	B	T
B 20/40	RED				B	R	B	L	L	R	B 20/40			T	T	L	B	T	R
C 20/30	AMBER				L	T	T	L	R	B	C 20/30			L	B	R	R	L	T
D 20/20	GREEN				R	L	T	L	B	T	D 20/20			R	T	R	L	T	B

**DIAL AT 10: COLOR PERCEPTION - FAR**    1-R    2-L    3-B    4-T    5-B    6-L    7-T    8-R

**DIAL AT 11: MUSCLE BALANCE - FAR**    Ball **IN** Box = **PASS**    Ball **OUT** of Box = **FAIL**

**DIAL AT 12: ALLEN PRE-SCHOOL FAR**

LINE/ACUITY	COLOR	BOTH EYES											
20/100		CAR, PHONE, CAKE, HORSE											
		LEFT EYE						RIGHT EYE					
A 20/50	BLUE	HORSE, CAKE, PHONE CAR						CAKE, HORSE, CAR, PHONE					
B 20/40	RED	CAKE, CAR, HORSE, PHONE						PHONE, CAKE, HORSE, CAR					
C 20/30	YELLOW	PHONE, HORSE, CAR, CAKE						HORSE, CAR, PHONE, CAKE					



1.800.344.9500 / 1.773.867.0380 / Fax: 1.773.867.0388  
 sales@stereooptical.com / www.stereooptical.com  
 8600 W. Catalpa Ave, Suite 703  
 Chicago, IL 60656-1116  
 © Stereo Optical Co., Inc. 2017